

<b>IHE Integration Statement</b>		<b>Date:</b>	November 2015
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	
Carestream	Touch Prime Ultrasound	1.0	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	Acquisition Modality	None	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Portable Data for Imaging (PDI)	Portable Media Creator (PMC)	None	
<b>Internet address for Carestream Health IHE information:</b> <a href="http://www.carestream.com/ihe.html">http://www.carestream.com/ihe.html</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>DICOM</b>	<a href="http://www.carestream.com/dicom.html">http://www.carestream.com/dicom.html</a>		
<b>Links to general information on IHE</b>			
North America: <a href="http://www.rsna.org/ihe.aspx">http://www.rsna.org/ihe.aspx</a>	Europe: <a href="http://www.ihe-europe.net/">http://www.ihe-europe.net/</a>	Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>	