

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Carestream health	Carestream PACS: Client Suite, Workflow Manager, IMS, CD Direct Suite, XDS Registry.	CSP v11.3	24 May 2011
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	Image Manager/Image Archive Image Display PPS Manager	None	
Access to Radiology Information (ARI)	Image Manager/Image Archive Image Display Report Repository Ext. Report Repository Access	None	
Patient Information Reconciliation (PIR)	Image Manager/Image Archive PPS Manager	None	
Consistent Presentation of Images (CPI)	Image Manager/Image Archive Image Display	None	
Cross-enterprise Document Reliable Interchange (XDR)	Document Source	None	
Cross-enterprise Document Sharing for Imaging (XDS.b-1)	Imaging Document Source	Set of DICOM Instances, Multiple documents None	
Cross-enterprise Document Sharing (XDS.b)	Document Repository Document Registry	Multiple document	
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer	None	
Evidence Documents (ED)	Image Manager/Image Archive	None	
Simple Image and Numeric Report (SINR)	Report Repository Ext. Report Repository Access	None	
Mammography Image (MAMMO)	Image Manager/Image Archive Image Display	None	
Key Image Note (KIN)	Image Manager/Image Archive	None	
Nuclear Medicine Imaging (NMI)	Image Manager/Image Archive	None	
Cardiac	Image Manager/Image Archive	None	

Catheterization Workflow(CCW)		
Audit Trail and Node Authentication (ATNA)	Secured Application Image Manager/Image Archive Report Repository Ext. Report Repository Access Image Display XDS-I Document source	None
Consistent Time (CT)	Time Client	None
Portable Data for Imaging (PDI)	Portable Media Creator Portable Media Importer Image Display Display	Web Content None None None
<b>Link to Carestream Health IHE information:</b> <a href="http://ihe.carestreamhealth.com">ihe.carestreamhealth.com</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>DICOM</b>	<a href="http://dicom.carestreamhealth.com">dicom.carestreamhealth.com</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>