

D-Speed



## Excellent images with fine detail and high contrast

Dental intraoral D-Speed film is a cost-effective, high-quality solution if you need the technique settings of a traditional dental radiographic system. Along with fine detail and high contrast, this film delivers outstanding processing and handling characteristics.

- A proven industry standard D-Speed film from Carestream Dental – a leader in the dental industry
- Processes in automated or manual equipment; no change in your current processing procedures or equipment is required
- Produces consistently high-quality images, ideal for diagnostic purposes
- Enhanced processing robustness – providing more tolerance to chemistry variations
- Thinner packets mean enhanced patient comfort
- Available in 100-packet boxes which allows you to use entire stock prior to expiration date
- Available in Size 2 to accommodate most dental exams

# Processing Information

For optimum results, process your radiographs using Carestream Dental chemicals.

## Recommended exposure settings\*

Settings		kV	50	60	65	65	65	70
		mA	7	7	7,5	8	10	7
Maxillary	Anterior		1.37	0.55	0.32	0.30	0.27	0.24
	Premolar		1.83	0.73	0.43	0.40	0.37	0.32
	Molar		2.06	0.82	0.48	0.45	0.41	0.38
Mandibular	Anterior		1.14	0.46	0.27	0.25	0.23	0.20
	Premolar		1.26	0.50	0.29	0.28	0.25	0.22
	Molar		1.37	0.55	0.32	0.30	0.27	0.24
Bite Wing	Premolar		1.37	0.55	0.32	0.30	0.27	0.24
Maxillary	Anterior		5.48	2.2	1.28	1.2	1.08	0.96
	Premolar		7.32	2.92	1.72	1.6	1.48	1.28
	Molar		8.24	3.28	1.92	1.8	1.64	1.52
Mandibular	Anterior		4.56	1.84	1.08	1	0.92	0.8
	Premolar		5.04	2	1.16	1.12	1	0.88
	Molar		5.48	2.2	1.28	1.2	1.08	0.96
Bite Wing	Premolar		5.48	2.2	1.28	1.2	1.08	0.96

\* NOTE: For large patients, increase time by approximately 25%. For children, small, or edentulous patients, decrease time by approximately 30%.

## Ordering Information for D-Speed

Size	Nominal Film Size	Invoice Description	Product Description	Films per Box	REF No
Size 2	30.5x40.5mm	D-Speed 100PKT 30.5x40.5 film	Vinyl packet, single film	100	105 1267

Detail is everything. That's why, despite having set the industry standard, we continue to search for ways to improve the image quality that our systems can deliver. It's all in the name of diagnostic excellence, which is fundamental to gaining your patients' confidence and optimizing their care.

For more information, contact your local authorized representative for Carestream Dental products or visit: [carestream.com](http://carestream.com).