

<b>IHE Integration Statement</b>		<b>Date:</b> .....July 2013
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>
Carestream Health Inc. ....	<i>DryView 7; 70 Laser Imager</i>	1.0
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
<b>Integration Profiles implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Consistent Presentation of Images (CPI)	Print Server	None
Mammography Image (MAMMO)	Print Server	None
<b>Internet address for Carestream Health IHE information:</b> <a href="http://www.carestream.com/ihe.html">http://www.carestream.com/ihe.html</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>DICOM</b>	<a href="http://www.carestream.com/dicom.html">http://www.carestream.com/dicom.html</a>	
<b>Links to general information on IHE</b>		
North America: <a href="http://www.rsna.org/ihe.aspx">http://www.rsna.org/ihe.aspx</a>	Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>

End of Document